

520 A 1500 – 400



Medical

- 520 A 1500 – 400 Ultimate
 - Part pump housing
 - Cavities 2
 - Cycle time 40s
 - Material LSR, Momentive 4070
 - Shot- / part weight 4g / 2g

520 A 1500 – 400

Highlight

- Application expertise LSR
 - Slider on middle platen
 - Sealing for vacuum around the sliders
 - Thermal separation cold channel to „hot“ mould
 - New Elmet dosing unit Top5100



520 A 1500 – 400

Highlight

- „Connected“
 - LSR dosing unit connected via EM 82.3
- „Automated“
 - Part removal and set down at passive cooling station with Yaskawa and EOAT from Barth



520 A 1500 – 400

- 520 A 1500 – 400 Ultimate

- Partner

- Customer

- PSG Dover, US

- Mold

- ELMET Elastomere Produktions- und Dienstleistungs GmbH, AT

- Dosing unit

- ELMET Elastomere Produktions- und Dienstleistungs GmbH, AT

- EOAT

- Barth Mechanik GmbH, D

- Material

- Momentive Performance Materials Company GmbH, D